

UNIVERSAL CR TECHNIQUE CHART 100% More mAs (225 Speed)

Part	View	Small		Medium		Large	
All techniques are male patients. If female, subtract 10 lbs.		120-160 lbs.		160-200 lbs.		200-240 lbs.	
		kV	mAs	kV	mAs	kV	mAs
Knee	AP/Lat (grid)	85	7.5	85	9	85	10.5
Knee	AP/Lat (non grid)	70	6	70	7.5	70	9
Knee	Sunrise (non grid)	70	7.5	70	9	70	10.5
L-Spine	AP (grid)	90	21	90	30	90	42
L-Spine	Oblique (grid)	90	33	90	48	90	66
L-Spine	X-Table Lat (grid)	95	48	95	66	95	90
Mandible	PA (grid)	85	9	85	12	85	15
Mandible	Obl (grid - 40")	81	6.9	81	9.6	81	12.3
Mandible	Obl (non grid - 40")	70	6	70	8.4	70	10.8
Nasal Bones	Lateral (non grid)	70	4.5	70	5.7	70	6.9
Pelvis	AP (grid)	85	18	85	27	85	36
Ribs	PA Upper (72")	85	18	85	24	85	30
Ribs	Obl Upper (72")	85	27	85	36	85	45
Ribs	AP Lower (40")	85	18	85	27	85	36
Shoulder	AP (grid)	85	9.3	85	13.5	85	18
Shoulder	AP (non grid)	70	6	70	8.1	70	9.9
Shoulder	Y View (grid)	85	25	85	34	85	48
Shoulder	Axillary (non grid)	70	6	70	8.4	70	10.5
Sinus	Caldwell	85	10.5	85	13.2	85	16
Sinus	Waters	85	12	85	15	85	18
Sinus	Lateral	85	5.4	85	7	85	8.7
Skull	PA	85	10	85	12.6	85	15.3
Skull	Lat (grid)	85	7.5	85	9	85	10.5
Tib-Fib	AP/Lat (non grid)	70	4.2	70	5.5	70	7.1
Toe	All Views	66	1.5	66	1.8	66	2.1
T-Spine	AP	85	16.5	85	22.5	85	28.5
T-Spine	Lat (breathing)	90	39	90	54	90	69
Wrist	PA	66	1.8	66	2.4	66	3
Wrist	Lat	70	2.4	70	3	70	3.6
Zygomatic Arch	AP Axial - Townes	80	7.8	80	9.6	80	10.8
Zygomatic Arch	SMV (grid)	80	7.5	80	9.3	80	11.1
Zygomatic Arch	SMV (non grid)	70	6.3	70	8	70	10